transcom

TRANSCENDENCIAS COMERCIALES, S. L.

C.I.F. B-48833404

TRANSCOM SL

COLON HYDROTHERAPY

HC – HEATER

OPERATING INSTRUCTIONS

transcom

TRANSCENDENCIAS COMERCIALES, S. L.

C.I.F. B-48833404

CONTENTS:

FUNCTION.....	
INSTALLATION	
CONNECTION TO MAINS WATER.....	
CONNECTION TO POWER SUPPLY	
PUTTING INTO OPERATION	
OPERATION OF THE DEVICE	
TEMPERATURE CONTROL	
DEVICE MAINTENANCE	
PREINSTALLATION DESCRIPTION STEP BY STEP...	
ILLUSTRATED SCHEME OF HC-HEATER	

transcom

TRANSCENDENCIAS COMERCIALES, S. L.

C.I.F. B-48833404

1. FUNCTION

The HC – HEATER device, is the newest machine serving as a connection for all type of intestinal rising machines known in the market.

The device incorporates a intermediate tank to have a mechanical separation between the water main line and the intestinal rising device. This system correspond to standard EN 1717.

2. INSTALLATION

The device has been designed with reduced dimensions, ideal to place under each standard patient bed.

The HC – HEATER should be place as close as possible to the colon device to prevent lost temperature in the transportation of the water.

3. CONNECTION TO MAINS WATER

First, the main water line must be connected. To connect it the device provides a connection kit, to make it easily.

In the same kit there are also, some spare parts to make the connection with the colon hydrotherapy device.

transcom

TRANSCENDENCIAS COMERCIALES, S. L.

C.I.F. B-48833404

4. CONNECTION TO POWER SUPPLY

The device works with 220/50Hz or 125/60Hz. Consumes 1300 w. The device goes with a spare parts kit to make the installation well and with some fuses to change when its owns broke.

The device is equipped with two switches to command. The first one is the general switch, this is the power switch of the device and the other switch is to heat the water when the deposit is full. Near the heater switch, the is a heating lamp. When the lamp turns off, means the water is hot and ready to work.

5. PUTTING INTO OPERATION

First, the device must be empty of water. This is important to carry it easily, without weight, and water.

To put it in operation, first the heater must be full of water. To make it better, press only the general switch, after installing all the connections. When the general switch is in on, the pump stars working filling the heater of water. This filling takes about 5 minutes. When the heater is full, so, when the water is going out the red tube, seems that the heater is ready to work.

We recommend to cover the cold water outlet when you are filling the heater and the red tube placing to a evacuate tube.

When the heater is full of water is ready to start working.

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