



transcom

TRANSCENDENCIAS COMERCIALES S.L.

*Colon hydrotherapy,
the therapeutic purpose of water*



We are a Spanish company that manufactures colon hydrotherapy equipment and disposable cannulas. Since 1980, when TRANSCOM SL was established, we have been able to evolve and remain in the vanguard of technology, which has led us to be market leaders and pioneers in the design and manufacture of automatic and electronic colon hydrotherapy equipment.

Via our extensive commercial network, we are present across all continents, and we have achieved a name and reputation at national and international level, as a result of which clinics, spas and therapy centers turn to us for projects and installations in their country.

With a continuous effort in R&D, we assure that our equipment is in line with the latest technology, which is justified by the quality approvals issued by the corresponding Ministries of Health.

*Colon Hydrotherapy,
the therapeutic
purpose of water*



Index

Company presentation 2



Index..... 3



Colon Hydrotherapy benefits 4



What is colon hydrotherapy? 5



Regeneration of the intestinal flora..... 5



HC-1 6



HC-2000 8



HC-3000 10



HC-1 "Classic" 12



HC-2000 "Classic" 13



Adult hydrokit 14



Small hydrokit 15



Straight hydrokit 16



Speculum + obturator..... 17



Gravity hydrokits..... 18



Hydrokit with limiter..... 19



Accessories 20



Colon Hydrotherapy

Benefits



Colon hydrotherapy helps us to prevent and improve all illnesses, and it is mainly indicated in the following pathologies:

- Constipation, diarrhoea
- Pregnancy (up to four months), beneficial for pregnant women
- Mucous, intestinal sluggishness, bloating, wind
- Oedemas, localized inflammatory conditions, ankylosis
- Dorsal cord pain, chronic lumbagos
- Intestinal parasitosis, candidiasis
- Dry skin
- Diverticulosis, colitis, Crohn's Disease
- Dermatitis, acne, eczema, psoriasis
- Depression, migraine, sleeping disorders, stress
- Insufficient renal elimination
- Cellulite, obesity, inadequate weight regulation
- Prevention of colon cancer
- Asthma, chronic bronchitis, sinusitis, rhinitis
- Weak abdominal muscle tone
- Vein and lymphatic disorders, tired and painful legs
- Cystitis, inflammations of the ovaries, dysmenorrhoea
- Stimulation of the body's natural defences
- Post-operative complications (after a general anaesthetic)
- Endocrinal disorders
- Sports medicine
- Recovery from strains
- Arthritis
- Hospital applications: Before performing an opaque enema or a rectoscopy and before performing surgery on the intestines
- Intestinal reduction in paraplegics and quadriplegics, ailing persons and in the presence of degenerative pathologies, all due to a chronic intoxication of the organism

What is colon hydrotherapy?

The cleaning of the intestines has been present since olden times. The first historical data can be found in Egypt, where the doctors-priests performed enemas or clysters for therapeutic purposes. After this, the Romans adopted this tradition, which has been maintained to the present day thanks to the appearance of technology specially designed to facilitate and improve the cleaning of the intestines.

Colon hydrotherapy is a therapy that has the objective of completely cleansing the colon, using water that has been previously filtered and purified, at a controlled temperature and pressure.

Since olden times, treatments have been performed using enemas of the large intestine. Colon hydrotherapy is a well-known possibility, which stems from the USA, where for many years its effectiveness has been accredited for cleansing and treating the intestine.

This process has become necessary as, nowadays, virtually no-one has an intestinal function in perfect conditions, and most people suffer from an altered composition of intestinal bacteria, which is known as dysbiosis.

As a result of this, the digestive processes are deficient leading to the formation of toxic substances and waste. In a new process the self-poisoning or autointoxication of the body is produced. The consequences of this are a lack of vitality, tiredness, depressions, lack of concentration, aggressiveness and states of anxiety. Infections, rheumatism, polyarthritis, hypertonia, skin conditions, migraines, allergies and other ailments are illnesses that are currently attributed to the incorrect functioning of the intestine.

Regeneration of the intestinal flora

Until the 19th century, oligoelements were considered to be organic impurities or analytical errors. Nowadays they are considered to be essential food complements for the functional integrity of human beings.

Twelve elements are assumed to be essential for the appropriate functioning of the human body: calcium, magnesium, potassium, sodium, zinc, cobalt, copper, chrome, iron, manganese, molybdenum and selenium.

The supply of oligoelements is carried out according to certain proportions between each one and at very low concentrations, as the more diluted the element the greater the contact surface of the complex with the organism. Their doses do not have the contraindications or side-effects of the same metals, metalloids and minerals at higher doses, as presented in some commonly used pharmaceutical specialities.

In short, given the circumstances that daily activity imposes on us, it is sensible to use the knowledge that science provides us in order not only to prevent specific illness but also to maintain the best state of health and wellbeing possible, and one way of achieving this could be to combine traditional pharmacopoeia with oligotherapy.





transcom



HC-1

Height: 1020/1100 mm.
Width: 550 mm.
Depth: 550 mm.
Weight: 35 kg.

*Useful and
easy-to-use*

**Essential for the correct cleaning
of the intestines.**


**The most highly-perfected manual
device in Europe.**

**With an efficient and safe
cleaning and disinfection system.**

**A practical, easy-to-use
equipment.**

Technical specifications:

- Innovated design
- Easy handling
- Easy installation
- Low energy
- Selector of voltage 210-220v
- Automatic pressure safety at 150 mbar
- Automatic temperature safety from 22° C to 40° C
- Digital thermostat
- Flow regulator from 0 to 100 l/h
- System of hermetic drainage not causing a foul smell
- Integrated deposit for disinfectant
- Polypropylene filters 5 μ c
- Introduction of therapeutic solutions
- UV lamp, sterilization of water
- Cleaning and disinfection system
- Regulable pressure gauge from 0 to 250 mbar
- Chromotherapy
- Height level regulable
- Incorporated wheels





transcom



HC-2000

Height: 1020/1100 mm.
Width: 550 mm.
Depth: 550 mm.
Weight: 41 kg.

*The benefits
of hydrotherapy
in an automatic
device*



**This device is the combination of
a refined, robust modern
technique, which is very easy to
use.**

**The complete cleaning of the large
intestine and unprecedented
safety are the basis of COLON
HYDROTHERAPY.**

Technical specifications

- Innovated design
- Easy handling
- Easy installation
- Low energy
- Selector of voltage 210-220v
- Automatic pressure safety at 150 mbar
- Automatic temperature safety from 22° C to 40° C
- Digital thermostat
- Flow regulator from 0 to 100 l/h
- System of hermetic drainage not causing a foul smell
- Integrated deposit for disinfectant
- Polypropylene filters 5 μ c
- Introduction of therapeutic solutions
- UV lamp, sterilization of water
- Cleaning and disinfection system
- Regulable pressure gauge from 0 to 250 mbar
- Irrigation liter counter, with lighted display
- Irrigation time counter, with lighted display
- A system of collection of samples for analysis
- Remote control
- Cromotherapy
- Height level regulable
- Incorporated wheels





transcom

HC-3000

Height: 1570 mm.
Width: 580 mm.
Depth: 580 mm.
Weight: 70 kg.


*High technology
at the therapist's
service*




Thanks to digital technology,
the HC-3000 incorporates
unprecedented electronic
safety and comfort systems,
guaranteeing the perfect
development of the Colon
Hydrotherapy.

Technical specifications

- Innovated design
- Easy handling
- Easy installation
- Low energy
- Automatic pressure safety at 150 mbar
- Automatic temperature safety from 22° C to 40°C
- Digital thermostat
- Electronic regulator of temperature
- Contrasts
- Flow regulator from 0 to 100 l/h
- System of hermetic drainage not causing a foul smell
- Integrated deposit for disinfectant
- Introduction of therapeutic solutions
- Polypropylene filters 5 μ c
- Automatic system for cleaning and disinfection
- Automatic pressure gauge from 0 to 250 mbar
- Irrigation liter counter, Irrigation time clock
- A system of collection of samples for analysis
- Incorporated ozone
- Wireless remote control
- Total shutdown system
- Incorporated wheels





HC-1 Classic


Height: 550 mm.
Width: 600 mm.
Depth: 210 mm.
Weight: 21 kg.

Technical specifications:

- Easy handling
- Easy installation
- Low energy
- Selector of voltaje 210-220v
- Automatic pressure safety at 100 mbar
- Automatic temperature safety from 22° C to 40° C
- Digital thermostat
- Flow regulator from 0 to 100 l/h
- System of hermetic drainage not causing a foul smell
- Integrated deposit for disinfectant
- Polypropylene filters 5µc
- Introduction of therapeutic solutions
- UV lamp, sterilization of water
- Cleaning and disinfection system
- Regulable pressure gauge from 0 to 250 mbar

HC-2000 Classic

Height: 750 mm.
Width: 600 mm.
Depth: 210 mm.
Weight: 29 kg.



Technical specifications:

- Easy handling
- Easy installation
- Low energy
- Selector of voltage 210-220v
- Automatic pressure safety at 100 mbar
- Automatic temperature safety from 22° C to 40° C
- Digital thermostat
- Flow regulator from 0 to 100 l/h
- System of hermetic drainage not causing a foul smell
- Disinfectant deposit
- Polypropylene filters 5 μ c
- Introduction of therapeutic solutions
- UV lamp, sterilization of water
- Cleaning and disinfection system
- Regulable pressure gauge from 0 to 250 mbar
- Irrigation liter counter, with lighted display
- Irrigation time counter, with lighted display
- A system of collection of samples for analysis
- Remote control




*High quality.
Maximum comfort*



- Material traslucid
Polypropylene
- Large size speculum
- Obturator perfect design
for an easier handling
- Waste helicoidal hose
manufactured in EVA
material 1300 mm
- Inflow hose 1500 mm
- Lubricant 2,5 gr of
vaseline

Small hidrokit



One-use kits



- Material traslucid
Polypropylene
- Small size speculum
- Obturator perfect design
for an easier handling
- Waste helicoidal hose
manufactured in EVA
material 1300 mm
- Inflow hose 1500 mm
- Lubricant 2,5 gr of
vaseline



Straight hydrokit

- Material traslucid
Polypropylene
- Straight speculum
- Obturator perfect design
for an easier handling
- Waste helicoidal hose
manufactured in EVA
material 1300 mm
- Inflow hose 1500 mm
- Lubricant 2,5 gr of
vaseline



Speculum + obturator



- Large size speculum
- Obturator perfect designed for an easier handling





Gravity hydrokits



- Material traslucid
Polypropylene
- Large size speculum
- Obturator perfect design
for an easier handling
- Waste helicoidal hose
with straight zones
- Inflow hose 1300 mm
- Lubricant 2,5 gr of
vaseline



Hydrokit with limiter



- **Material traslucid Polypropylene**
- **Speculum with limiter**
- **Obturator perfect design for an easier handling**
- **Waste helicoidal hose manufactured in EVA material 1300 mm**
- **Inflow hose: PVC 1500 mm long (optional 6 mm / 8 mm / 10mm diameter)**
- **Lubricant 2,5 gr of vaseline**





Disinfectant

In order to complete a Colon Hydrotherapy session you can find a whole range of complements and accessories:

- Bottle to insert therapeutic solutions. - Hydro-colon couch: Specially designed for colon hydrotherapies.
- Mini-kits: Disposable catheter for the MH-1 irrigator or enema device.
- Disinfectant: For quick, simple and safe high-level disinfection (HDL), maintaining excellent compatibility with all materials.
- Poster: Diagrams of the colon.

Besides all the products listed above, you can find a full range of disposable goods to be used in one session, such as:

Disinfectant

Filters HC-1

Filters HC-2000

Filters HC-3000

Under pads

Gowns

Gloves

Disposable slippers

Sheets and Hydro-colon couch



Filters HC-1 / HC-2000

Filters HC-3000





Disposable slippers



Under pads

Oligoelementos de origen natural



Gowns




Gloves



*Complements and accessories
in order to complete
a Colon
Hydrotherapy session*

Sheets and Hydro-colon couch





transcom

Distributors

Bulgaria
China
France
Indonesia
Iran
Italy
Japan
Jordan
Kazakhstan
Kuwait
Lithuania
Mexico
Portugal
Moldova
Romania
Russia
Spain
South Africa
Syria
Thailand
Turkey
Ukraine
United Arab Emirates
United Kingdom
United States of America
Vietnam





"The disease starts and the health begins in the colon"



FACTORY | TRANSCOM

Parque empresarial Zuatzu - Edificio Urumea
C/ Zubiberri, 31 - 20018 San Sebastián - Gipuzkoa (Spain)
www.transcomsl.com - e-mail: info@transcomsl.com
Tel. (34) 943 224 360 - Fax (34) 943 224 275